

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS20 4370

GLDP | Glycine decarboxylase P protein

Product information

Immunogen | KLH-conjugated peptide, highly conserved in Arabidopsis thaliana GLDP1 UniProt: Q94B78, TAIR: At4g33010 and

GLDP2 UniProt: <u>O80988</u>, TAIR: <u>At2g26080</u>

Host Rabbit

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 50 ul

Reconstitution For reconstitution add 50 μl, of sterile water

Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please Storage remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

Application information

Recommended dilution 1:50 (IL), 1:5000 (WB)

Expected | apparent MW

113 kDa

Confirmed reactivity Arabidopsis thaliana, Flaveria spp., Neurachne spp., Oryza sativa, Zea mays

Predicted reactivity grasses, non-succulent dicots

Species of your interest not listed? Contact us


Not reactive in Tecticornia spp., Salsola soda

Additional information So far only leaf tissue has been used for immunolocalization and western blot experiments

Khoshravesh et al. (2020). The Evolutionary Origin of C4 Photosynthesis in the Grass Subtribe Neurachninae. Plant Selected references

Physiol. 2020 Jan;182(1):566-583. doi: 10.1104/pp.19.00925.

application example



6 µg of total protein extracted freshly from Z. mays = Zea mays, A. tha= Arabidopsis thaliana, O. sat = Oryza sativa, was freshly extracted from leaf tissue in extraction buffer (50 mM HEPES-KOH pH 7.1, 1 mM EDTA pH 8.0, 10 mM DTT, 5 mM MgCl2, 1 mM PMSF, 20 μg/ml chymostatin, 25 µg/ml aprotinin, 1X Protease Inhibitor Cocktail VI, Plant Cell (AG Scientific)). Proteins were denatured with sample buffer (1X = 62.5 mM Tris-HCl pH 6.8, 2% [v/v] SDS, 10% [v/v] glycerol, 0.01% [w/v] bromophenol blue (Sigma), 100 mM DTT), heated to 95°C in a water bath for 3 minutes and collected at 13523 x g in a bench centrifuge. 6µg of protein was loaded per well and separated on a 12% resolving gel (0.375 M Tris-HCl pH 8.8, 12% [v/v] acrylamide, 0.1% [v/v] ammonium persulfate and 0.05% [v/v] TEMED) and 5% stacking gel (0.125 M Tris-HCl pH 6.8, 5% [v/v] acrylamide, 0.1% [v/v] SDS, 0.1% [v/v] ammonium persulfate and 0.025% [v/v] TEMED). Proteins were blotted onto nitrocellulose membrane (pore size 0.45 μm) using Trans-Blot Semi-Dry Transfer Cell (Bio-Rad), and transfer was run at 25 mA per gel stack for 1½ hours. Blots were blocked in blocking solution (5% [w/v] skim milk powder in TBST (12.5 mM Tris, 137 mM NaCl, 2.7 mM KCl, 0.1% [v/v] Tween-20 detergent)) overnight at room temperature with gentle agitation. The blot was rinsed three times quickly in TBST, then incubated in the primary antibody solution at a dilution of 1:5000 in blocking solution for 1 hour at room temperature with gentle agitation. Blots were rinsed three times quickly in TBST, then washed three times for 5 minutes in TBST at room temperature with gentle agitation. Blots were then incubated in secondary antibody solution (goat anti-rabbit IgG horse radish peroxidase (Sigma) diluted 1:3000 in blocking solution) for 1 hour at room temperature with gentle agitation. The blot was washed three times for 5 minutes in TBST at room temperature with gentle agitation, then rinsed twice in TBS (12.5 mM Tris, 137 mM NaCl, 2.7 mM KCl). Blots were developed with chemiluminescent detection reagent, following manufacture's



This product is **for research use only** (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

recommendations, incubating at room temperature for 1 minute and imaged with Bio-Rad ChemiDoc system. Exposure time was 40 seconds.

Courtesy of Prof. Martha Ludwig, School of Molecular Sciences, Australia